



Save up to 40 % with the anniversary offers
on innovative MELAG solutions

Follow us :    

www.melag.com/en

System solution, innovation and quality

Made in Germany. Made for You.

Would you like to save time and money on instrument reprocessing? With our system solutions, we ensure efficient and safe workflows of the highest quality. As an owner-operated family company founded in 1951, we specialize in products aiming at optimized hygiene in practice and clinics environments. Get an overview of our advantages:



Dr. Niklas Gebauer and Sebastian Gebauer, CEOs



Experience
Specialization since 1951



Quality
100 % made in Germany



Innovation
500 specialized employees in Berlin



Trust
World-leading with over 700,000 satisfied customers



System solutions
Instrument reprocessing aligned to your specialism



In this video, we introduce you to the bright minds behind the Careclave-project. Dedicated employees who always have the health of people around the world in focus.

You have questions?
We the answers!

Pete Gibbons

Mobile : +44 7961 439355

E-Mail: gibbons@melag.eu



MELAcare Deal

NEW



The game-changer for handpiece cleaning, disinfection and lubrication: MELAcare 600 02

MELAcare & Vacuclave 105 Deal



The dream team for ultra-fast instrument reprocessing: MELAcare 600 & Vacuclave 105 04

Vacuclave 118 Deal



The MELAG price-performance champion: Vacuclave 118 06



MELAcare

Wellness for your handpieces and turbines

Do you know that feeling when a new dental handpiece just works perfectly? When treatment goes smoothly and everything runs like clockwork? That feeling now has a name: MELAcare – for fully automatic cleaning, disinfection and oil care with many benefits:

- ✓ MELAcare Blue lid for cleaning, disinfection and oiling of up to 6 handpieces
- ✓ MELAcare Green Lid for cleaning and disinfecting of e.g. 6 ultrasonic scaler tips
- ✓ Reprocessing procedure in accordance with EN ISO 15883-1/2 from approx. 14 minutes
- ✓ Gentle on instruments and chemical-free thanks to ULTRAflow Technology
- ✓ Precise lubrication and chuck care for reliable functionality and a long service life of your handpieces
- ✓ Large intuitive touch screen display for ease of use
- ✓ Smart documentation and approval directly on the included USB stick
- ✓ Incredibly compact design with a small footprint



MELAcare bundle:

1x MELAcare 600,
1x MELAcare Lid Blue,
6x Adapter (individually configurable for your instruments)
1x MELAdem 40

Anniversary Offer : **40 %**
£2,990 +VAT

Vacuclave 105 & MELAcare

The fastest combo for handpiece and instrument reprocessing

Reprocessing handpieces is not only essential for patient safety – it is also a key economic factor in everyday practice. Fast workflows with reliable reprocessing and care results can reduce the need for additional instruments while helping to protect valuable handpieces from downtime, wear and costly repairs in the long term.

- ✓ **Vacuclave 105:** The compact 5-litre autoclave for sterilising up to 22 handpieces or up to 2kg of instruments from 6.5 min
- ✓ **MELAcare:** The compact unit for cleaning, disinfection and care of up to 6 handpieces from 14 min*
- ✓ **Long-term instrument protection** through brilliant cleaning results without chemicals, combined with precise oiling and chuck care
- ✓ **Ultra-fast reprocessing cycles** increase throughput, thereby reducing the need for additional costly instruments.



Vacuclave & MELAcare bundle:

- 1x Vacuclave 105
- 2x Tray flat
- 1x MELAcontrol Helix
- 1x MELAcare 600
- 1x MELAcare Lid Blue
- 6x Adapter (individually configurable for your instruments)

Anniversary Offer : **-38 %**
£4,990 +VAT

Speed up your practice workflow – from instrument care to patient safety

14 min for cleaning, disinfection and, where applicable, care
+ 6,5 min for strerilisation



* Cycle time with empty chamber and optimum conditions

Pro *Line* : Vacuclave 118

Discover the new MELAG price-performance champion!

You will love the all-new MELAG device generation - from the functions to the price! Because the Pro Line stays true to its line: Outstanding performance, smart functions as well as a modern and extremely robust design made in Germany:

- ✓ **Powerful Class B vacuum process** for up to 6 kg on 5 trays
- ✓ **Smart-Touch Display** for easy operation just like your smartphone
- ✓ **Power Save and Easy Care Concepts** for **sustainable and user-friendly reprocessing workflows**
- ✓ **Compliance with legal requirements** through logging via **USB or network**



Vacuclave 118 bundle:

- 1x Vacuclave 118, USB stick, LAN network connection
- 3x Trays
- 1x A Plus tray holder
- 1x MELAcontrol Helix
- 1x MELAdem 40 and water connection set
- 1x MELAJet

Anniversary Offer : **-25 %**
£3,749 +VAT

Notes on additional accessories:

[illegible]

Notes on additional accessories:

[illegible]

I hereby place an order:

For the bundles:

☐ MELAcare 600

☒ MELAcare Lid Blue with Adapters

☒ MELAdem 40

☐ Vacuclave 105 & MELAcare 600

☒ with 2 Trays

☒ MELAcontrol Helix

☒ MELAcare Lid Blue with Adapters

☐ Vacuclave 118

☒ Mount with 3 Trays

☒ MELAdem 40

☒ MELAjet

☒ MELAcontrol Helix

MELAG75

**You have questions?
We have the answers!**

Pete Gibbons

Mobile: +44 7961 439355

E-Mail: gibbons@melag.eu



Practice name: _____

Retail partner: _____

Practice address: _____

Date,
signature: _____

Offer is valid until 31.12.2026. All prices excl. 20 % VAT. Prices excl. costs for installation and commissioning. Joint action of participating partner and MELAG.